



# Thermometers for dilution fridges (V.Zavjalov, 11.2018)

The goal is to make standard thermometers for use in dilution fridges. I use carbon OHMITE resistors (OD/OF series) [1], which are claimed to be good thermometers at millikelvin temperatures [2]. The series OD and OF have different sizes (diameter 2.5 mm and 3.8 mm and length 7 mm and 10 mm respectively). I grind resistors down to approximately 0.7 mm thickness and glue it in a slit in a copper box using Stycast 2850. This provides a good thermal contact with the box, mechanical stability, and protection against moisture.




Copper boxes are ordered in Schaeffer-AG company [3]. The box has 30 × 30 × 5 mm size. It contains a 20 × 18 × 3 mm compartment for RC lowpass filters, a 1 × 16 × 4 mm slit for the resistor and a 16 × 7 × 3 mm additional compartment which can be used for sample chips. Holes are compatable with standard 20 × 20 M4 hole grid on BlueFors refrigerator plates. Drawings in FPD format are available [4], price was 315€ for 10 boxes and 10 lids.



PCBs for filters are ordered in Multi-CB company [5]. Price is 75€ for 10 boards with two different filters and two small sample plates each. Kicad project [6] and Gerber files [7] are available.

To connect thermometers to BlueFors 4-pin sockets we use inserts of LEMO FFA.0S.304.CLAC44 connectors (Farnell number: 2442870).



## List of devices

N	Resistor	$R_0$	$R_1$	$R_2$
2018-11-02 N1	OD270JE	27	56	56.0
2018-11-02 N2	OD270JE	27	59	64.4
2018-11-02 N3	OF121JE	120	231	277

Here  $R_0$ ,  $R_1$ , and  $R_2$  are nominal resistance, resistance after grinding, and resistance after soldering and glueing. Change during glueing probably means that contacts with carbon in the resistor after grinding are not stable.

## References

[1] Ohmite OD/OF series datasheet,  
[https://www.ohmite.com/assets/docs/res\\_od\\_of\\_oa.pdf](https://www.ohmite.com/assets/docs/res_od_of_oa.pdf)

[2] N. Samkharadze, A. Kumar, G. A. Cs  thy, A New Type of Carbon Resistance Thermometer with Excellent Thermal Contact at Millikelvin Temperatures, *JLTP*, **160**, 246–253 (2010),  
<https://doi.org/10.1007/s10909-010-0192-5>

[3] Schaeffer-AG company,  
<https://www.schaeffer-ag.de/en/>

[4] Box drawing in FrontPanel Designer format,  
[https://github.com/slazav/he3notes/raw/master/20181105-mk\\_therm/suppl/box\\_v1.zip](https://github.com/slazav/he3notes/raw/master/20181105-mk_therm/suppl/box_v1.zip)

[5] Multi-CB company,  
<https://portal.multi-circuit-boards.eu>

[6] Filter PCB, KiCAD project,  
[https://github.com/slazav/he3notes/tree/master/20181105-mk\\_therm/suppl](https://github.com/slazav/he3notes/tree/master/20181105-mk_therm/suppl)

[7] Filter PCB, Gerber files for ordering,  
[https://github.com/slazav/he3notes/raw/master/20181105-mk\\_therm/suppl/filter\\_pcb\\_v1.zip](https://github.com/slazav/he3notes/raw/master/20181105-mk_therm/suppl/filter_pcb_v1.zip)